


Referenzwerteliste Skyla Exoten (SI - Einheiten)



Parameter	Einheit	Amazonen- papagei blau	Amazonen- papagei gelb	Edel- papagei	Grau- papagei	Kakadu (grey cheek)	Molukken- kakadu	Weiß- hauben- kakadu	Wellen- sittich	Nymphen- sittich	Kanarien- vogel	Sittich	Ara Blau & Gold	Hyazinth- Ara	Hellroter Ara	Schlange	Schild- kröte	Meeres- schild- kröte	Eidechse
ALB	g/L	7-17	7-13	6-13	8-49	8-49	10-29	9-27	8-49	8-49	8-49	8-49	4-19	14-20	6-21	N/A	13-30	10-25	N/A
ALP	U/L	15-109	31-198	97-668	24-94	0-3000	40-160	30-204	54-326	0-3000	0-3000	0-3000	27-205	3-84	9-74	80-145	36-156	24-48	60-99
ALT	U/L	0-12	N/A	0-10	0-14	0-631	0-1	0-7	0-13	0-631	0-631	0-631	0-4	0-1	N/A	3-242	N/A	N/A	N/A
AMY	U/L	402-882	264-798	425-981	211-519	0-2500	414-972	424-759	187-582	0-2500	0-2500	0-2500	295-506	283-595	217-710	N/A	N/A	N/A	N/A
AST	U/L	73-160	67-151	19-159	28-200	150-400	53-140	84-143	55-154	100-350	150-350	125-350	17-126	68-122	32-95	5-35	14-18	66-315	5-103
BUN	mmol/l	N/A	N/A	N/A	1,1-41,8	1,1-41,8	N/A	N/A	1,1-41,8	1,1-41,8	1,1-41,8	1,1-41,8	0-1,1	0-1,4	0-2,5	0-3,2	6,4-12,1	12,5-37,1	0-3,6
Ca	mmol/l	2,03-2,98	1,93-2,75	1,98-3,15	1,7-2,43	0,48-4,1	1,83-3,23	2,1-2,9	1,53-2,95	2,15-3,48	0,48-4,1	2-4,08	1,75-3,35	1,85-2,63	2-2,8	2,58-6,15	2,58-3,93	1,1-3,78	1,88-2,58
CPK	U/L	97-450	98-646	144-418	71-408	100-300	60-606	112-559	54-252	100-300	0-2015	100-300	62-393	284-871	76-406	N/A	N/A	N/A	N/A
CREA	µmol/l	0-9	0-9	9-18	9-35	9-35	9-53	9-18	9-35	9-35	9-35	9-44	9-18	9-27	9-11	9-44	9-35	9-27	9-18
FRU	µmol/l																		
GGT	U/L	1-12	2-19	0-1	0-900	0-900	0-1	0-3	0-900	0-900	0-900	0-900	0-1	0-1	0-1	1,0-1,5	N/A	N/A	1-10
GLU	mmol/l	12,9-16,4	12,2-16,4	12,9-21,8	12,4-17,1	11,1-19,4	10,5-14,5	10,7-15,4	14,1-22,1	11,1-25	11,1-25	11,1-19,4	12,6-17,2	11,9-15,3	11,4-16,2	0,6-6	N/A	4-6,9	3-11
LAC	mmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
LDH	U/L																		
LIPA	U/L	0-25	0-25	0-25	0-35	0-35	0-25	0-25	0-35	0-35	0-35	0-35	0-30	0-30	0-30	N/A	N/A	N/A	N/A
Mg	mmol/l	1,32-1,56	0,95-1,52	0,74-1,28	0-1,56	0-1,58	0,7-1,44	0,74-1,28	0-1,56	0-1,58	0-1,56	0-1,56	0,82-1,56	0,53-1,15	0,62-1,32	N/A	N/A	2,76-3,87	N/A
PHOS	mmol/l	0,7-1,2	0,5-1	0,8-1,6	0,4-1,5	0,3-3,8	0,6-1,5	0,3-1,6	0,4-0,7	0,3-3,8	0,3-3,8	0,3-3,8	0,6-2	0,6-1,5	0,6-1,6	0,9-1,6	0,8-3,2	1,6-2,6	0,7-1,5
BA	µmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
TBIL	µmol/l	N/A	0-1,7	0-1,7	0-458,4	0-458,4	0-3,4	0-5,1	0-458,4	0-458,4	0-458,4	0-458,4	0-17,1	0-3,4	0-5,1	N/A	1,7-10,3	N/A	N/A
CHOL	mmol/l	5,4-11,3	3,8-7,8	4,8-8,3	5,6-8,1	0,7-11,7	3,6-8,3	3,2-9,0	4,5-7,1	0,7-11,7	0,7-11,7	0,7-11,7	2,5-5,0	2,2-4,2	3,1-7,9	1,8-3,8	N/A	4-10,7	1,7-3,8
TRIG	mmol/l	1,2-3	1-2,4	0,8-3,4	0,6-1,6	0-4,2	1,2-3,3	2,2-4,2	1,2-3,1	0-4,2	0-4,2	0-4,2	1-3,3	0,9-2,3	0,9-3	0,6-2	N/A	0-2,9	0,6-1,2
TP	g/L	37-50	34-46	35-50	26-49	25-45	36-55	38-55	20-30	30-50	30-50	25-45	37-52	39-46	38-51	29-80	30-70	36-68	30-81
T4	nmol/L	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
tCO2	mmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
UA	mg/dl	101,1-368,8	60-393	167-553	208-446	262-750	179-797	226-732	202-541	232-690	262-750	149-214	65-857	232-1231	125-506	83-631	95-190	42-131	184-506
#GLOB	g/dl																		
#UPC																			
#UREA	mmol/l																		
Na	mmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
K	mmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Cl	mmol/l	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

Referenzwerte Skyla Exoten Pro0474la, Vers. 1.0.xlsx