

SARS-COV-2 LYSIS KIT

RAPID RELEASE OF RNA FOR RT-PCR VIRAL ASSAYS



RNA release extraction-free protocol

The Leasy kit is designed to enable a fast and efficient release of viral RNA from human nasopharyngeal swabs in Universal Transport Medium™. The recovered RNA is suitable for TaqMan-based Real Time PCR assays.



Leasy Kit is intended for research use only,
not for use in diagnostics procedures

siliconbiosystems.com



MENARINI
silicon biosystems



RELEASE OF RNA STARTING FROM:

- human nasopharyngeal swab samples collected in Universal Transport Medium (UTM®).



FAST

UTM® to RT-PCR in less than 40 minutes



MINIMAL HANDS ON

just two pipetting steps



VERSATILE

compatible with TaqMan RT-PCR kits using exogenous internal control



HIGH THROUGHPUT

suitable for automated liquid handling platforms

Leasy kit is intended for research use only, not for use in diagnostics procedures.

Distributed by



Leasy kit is manufactured by



menariniagnostics.com

©2020 All rights reserved.

CMD_MKT_033 Rev1, Dec 2020